

Dale®

Hold-n-Place® Catheter Securement

Application Instructions

1. For IV and arterial catheters the Luer connection should be attached to the inserted catheter.
2. Clean & Prep the skin per hospital protocol.
3. On some patients hair removal or skin prep may be required for better adhesion.
4. Make sure the skin is clean and dry before applying the securement device.
5. Pull the top adhesive **PULL** tab back to the perforations on the two strips. **(A)**
6. Slide the securement device under the catheter/sheath positioning where the catheter tube joins the hub is centered over the extended point in the center of the base. **(B)**
7. Apply gentle pressure to the top of the hub to secure in the exposed adhesive. **(B)**
8. Hold the securement device and the catheter in position while pulling one bottom release liner and press adhesive to skin, repeat for other side. **(C)**
9. Remove the **PULL** tab completely to expose the adhesive strips. **(D)**
10. Secure one adhesive strip over and across the catheter hub applying gentle pressure to the top side of the adhesive strip, then repeat with the other adhesive strip to achieve a chevron securement. **(E)**
11. Press down on securement device to ensure better adhesion.
12. Apply a transparent dressing over the insertion site. You can place the transparent dressing over the entire securement device, however, only mid-way covering the securement device is recommended to ease dressing changes.
13. Follow hospital protocol for dressing change interval. If soiled or saturated in fluid, replace securement device.

Tape-Strips are included to help stabilize line or tube while applying the securement device.

Removal

Alcohol swab may be used to help with removal.

Ordering Information

848	Hold-n-Place IV, Arterial and Mid-line Catheter Securement with small Prevhex™ CHX transparent dressing	50/box
------------	---	--------

854	Hold-n-Place CVC, PICC and Arterial Sheath Securement with large PrevaHex™ CHX transparent dressing	50/box
------------	---	--------

For single patient use. Disposable.

For an application video please visit www.dalemed.com

Dale medical products are available from your preferred hospital/medical supply dealer.

Not made with natural rubber latex.

Patent www.dalemed.com/patents

Dale® & Hold-n-Place® are registered trademarks of Dale Medical Products, Inc. in the USA & EU.

PrevaHex™ is a trademark of entrotech life sciences Inc.

STERILE R

TEL.: +1-508-695-9316
www.dalemed.com

