

THINK

DALE®

Dale® Tracheostomy Tube Holders

240 Blue® & PediPrints®

EDUCATION
INNOVATION
PERFORMANCE
RELIABILITY
VALUE

The **Dale Tracheostomy Tube Holders** are designed to provide secure positioning and to minimize movement of the tracheostomy tube for all your patients. Whether you choose the well-known Dale® 240 Blue® for adults or PediPrints® series for kids, you can be sure it's the right choice.



Dale® 240 Blue®

- Adjustable neckband with exclusive moisture repellent lining
- Secure and easy to use fastener tabs fit any size trach plate
- Stretch material accommodates cough and edema
- 1" wide band fits up to 19.5" neck



Dale® PediPrints®

- Adjustable neckband with exclusive moisture repellent lining
- Printed neckband
- Secure and easy to use fastener tabs fit any size trach plate
- 241 - 1" wide band fits up to 18" neck
- 242 - 3/4" wide band fits up to 9" neck

THINK

DALE®

Dale® Tracheostomy Tube Holders

240 Blue® & PediPrints®

EDUCATION

INNOVATION

PERFORMANCE

RELIABILITY

VALUE

Unnecessary movement of the trach tube can result in accidental dislodgement or displacement of the tube, tracheopharyngeal fistula, tracheal stenosis, or airway granuloma. That's why it is important to know that no other securement method stabilizes tracheostomy tubes better than Dale®. Now you can eliminate the frustrations associated with twill ties and other holders while minimizing secondary complications at the same time.



Key Features

- Secure fastener tabs limit movement of the tracheostomy tube reducing accidental decannulation
- Exclusive moisture repellent neckband lining reduces the risk of skin breakdown by minimizing the amount of moisture absorbed
- Adjustable two piece neckband
- PediPrints® series offers kids fun, comforting prints†
- Fastener tabs quickly and easily attach to any size tracheostomy tube
- Stretch material‡ allows for cough reflex, accommodates edema, and ensures a snug fit without compromising security

Ordering Information

Product #	Description
#H84102401	Dale® 240 Blue® Tracheostomy Tube Holder 1" wide band fits up to 19.5" neck
#H84102411	Dale® 241 PediPrints® Tracheostomy Tube Holder 1" wide band fits up to 18" neck
#H84102421	Dale® 242 PediPrints® Tracheostomy Tube Holder 3/4" wide band fits up to 9" neck

Individually packaged. For single patient use.
Sold in boxes of 10

HSPCS Code: A7526

Latex-Free. In response to clinical concerns regarding patient hypersensitivity, Dale products contain no latex.

Dale products are available from your preferred hospital/medical supply dealer.

Use in a supervised care setting only. Pediatric, confused and other patients requiring special care should be monitored to prevent unintentional disengagement.

† On the 241 and 242 tube holder
‡ On the 240 Blue® tube holder

800-343-3980 • www.thinkdale.com

Dale®, 240 Blue®, and PediPrints are registered trademarks of Dale Medical Products, Inc. in the USA & EU.
©2011 Dale Medical Products, Inc.

Rev G ps-TTH

