

THINK

DALE®

Dale® ACE* Connector®

*(Access Controller for Enteral)

Closed system connector with needleless syringe port seal

Reduces Risk of Infection for You and Your Patient

100% closed-system Universal Connector provides in-line feeding, suctioning, irrigation, and medication delivery, reducing the risk of potentially infections splash back and cross contamination.



EDUCATION

INNOVATION

PERFORMANCE

RELIABILITY

VALUE

THINK

DALE®

Dale® ACE* Connector®

*(Access Controller for Enteral)

Closed system connector with needleless syringe port seal

EDUCATION

INNOVATION

PERFORMANCE

RELIABILITY

VALUE



The Dale ACE (Access Controller for Enteral) Connector protects you from potentially infectious splash back while shielding your patient from external contamination. This 100% closed system connector controls fluid flow into and out of Salem Sump®, Levin, and PEG/g-tubes, reducing the risk of cross contamination. The Dale ACE Connector replaces 5-in-1 connectors and enteral feeding Y-port adapters to provide in-line feeding, suctioning, irrigation, and medication delivery without the need to disconnect. The connector system remains closed to the atmosphere at all times for optimal infection protection.

Features and Benefits

- A convenient, compact device for controlling fluid flow into or out of a tube
- Connects to all Salem Sump®, Levin, NG tubes, other stomach tubes, and enteral feeding PEG tubes. One size fits all.
- Replaces enteral feeding tube "Y" adapters
- Replaces 5 in 1 connectors
- Protects the clinician from potentially infectious splash events by maintaining a closed system (no disconnect)
- Reduces risk of cross contamination, eliminates mess and formula loss by maintaining a closed system (no disconnect)
- Has a simple ON/OFF handle fluid controller
- Has a built in Needleless Syringe Port Seal



Misconnection Prevention Safety Feature

- Designed to physically prevent the administration of enteral formula into an I.V. line
- Incompatible with luer connectors
- Meets ANSL/AAMI ID 54:1996 standards

Ordering Information

Product #	Description
#H84104751	Dale® ACE* Connector *(Access Controller Enteral)

Individually packaged. For single patient use.
Sold in boxes of 50

US Patent #6,165,168

Latex-Free. In response to clinical concerns regarding patient hypersensitivity, Dale products contain no latex.

Dale products are available from your preferred hospital/medical supply dealer.

800-343-3980 • www.thinkdale.com

Dale® and ACE Connector® are registered trademarks of Dale Medical Products, Inc. in the USA & EU. Salem Sump® is a registered trademark of the Covidien AG Corporation in Switzerland. ©2011 Dale Medical Products, Inc.

Rev C ps-475

