

THINK

DALE®

Dale® Nasal Dressing Holder

The **Dale Nasal Dressing Holder** replaces adhesive tape to hold nasal dressings in place following rhinoplasty, septoplasty, sinus surgery, facial trauma, and nose bleeds.



Key Benefits Include:

- >> Totally adjustable
- >> Eliminates skin irritation
- >> Simplifies dressing changes

800-343-3980 • www.thinkdale.com

EDUCATION
INNOVATION

PERFORMANCE

RELIABILITY

VALUE

THINK

DALE®

Dale® Nasal Dressing Holder

The Dale Nasal Dressing Holder replaces adhesive tape to hold nasal dressings in place following rhinoplasty, septoplasty, sinus surgery, facial trauma and nose bleeds.

It consists of a clear plastic Dressing Holder, and is secured by ear loops. Its hook and loop fastener tabs can be infinitely adjusted to fit any patient comfortably. The soft material of the Dressing Holder adds to the comfort, while the plastic center section holds dressings in place without absorbing secretions. The Nasal Dressing Holder is easy to remove and replace, making dressing changes simple and painless. And because it requires no adhesive tape, the Nasal Dressing Holder won't damage skin.



Key Features

- Eliminates irritation and excoriation of facial skin
- Improves patient comfort
- Simplifies dressing changes
- One size fits most patients
- Saves nursing time
- Latex-Free

Ordering Information

Product #	Description
#H84106001	One size fits most Sold in boxes of 10

For single patient use.

Latex-Free. In response to clinical concerns regarding patient hypersensitivity, Dale products contain no latex.

Dale products are available from your preferred hospital/medical supply dealer.

EDUCATION
INNOVATION

PERFORMANCE

RELIABILITY

VALUE

800-343-3980 • www.thinkdale.com

Dale® is a registered trademark of Dale Medical Products, Inc. in the USA & EU.
©2011 Dale Medical Products, Inc.

Rev B ps-600

