



## Tracheostomy Tube Holder

### *CLINICAL ADVANTAGES*

<b>Features</b>	<b>Benefits</b>	<b>Clinical Advantages</b>
<b>Fastener Tabs</b> ·Fungus Resistant ·Tapered	<ul style="list-style-type: none"> <li>• Ensures hygiene and cleanliness</li> <li>• Helps prevent movement of the trach plate</li> <li>• Fits all brands of tracheostomy tubes with flanged ends</li> <li>• Narrow fastener tab easily threads through trach plate</li> <li>• Fastener tabs can be fastened back onto self to minimize absorption of secretions on a fresh trach or can be fastened back onto the neckband.</li> <li>• Securely holds trach tube in place</li> <li>• Will hold even if tabs are wet</li> </ul>	<ul style="list-style-type: none"> <li>• Increases caregivers' confidence and sense of security</li> <li>• Fungus resistant barrier reduces likelihood of nosocomial infection</li> <li>• Helps to keep stoma area clean and dry</li> <li>• Limiting movement at trach plate reduces likelihood of skin breakdown</li> <li>• Limiting movement of tracheostomy tube limits tracheal irritation by the tube, preventing damage to the tracheal mucosa</li> <li>• Reduces the risk of accidental decannulation and need for reintubation</li> <li>• Enhances recovery by assuring intact airway</li> <li>• Allows for easy change of tube holder and trach dressings</li> <li>• Prevents complications that can increase length of hospital stay</li> </ul>
<b>Moisture Resistant Neckband</b>	<ul style="list-style-type: none"> <li>• Keeps area dry</li> <li>• More comfortable than foam or twill ties</li> <li>• Generous one inch width prevents rolling or folding</li> <li>• Gentle, even on delicate skin</li> </ul>	<ul style="list-style-type: none"> <li>• Prevents skin excoriation on neck</li> </ul>
<b>Adjustable Neckband Closure</b>	<ul style="list-style-type: none"> <li>• More secure than twill ties</li> <li>• Easier and faster to change than twill ties</li> <li>• No cutting required means no fear of accidentally cutting neck or tube's inflation stem</li> <li>• Reduces staff time required for tracheostomy tie replacement</li> <li>• Adjusts to fit most patients, from pediatric to adult</li> </ul>	<ul style="list-style-type: none"> <li>• Promotes better skin care and prevents chaffing and skin breakdown</li> <li>• Allows quick and easy adjustment of neckband</li> </ul>
<b>Stretch Material on Neckband</b>	<ul style="list-style-type: none"> <li>• Accommodates weight of ventilator circuitry</li> <li>• Maintains a snug, comfortable fit</li> <li>• Accommodates edema and cough reflex</li> </ul>	<ul style="list-style-type: none"> <li>• Will not bind or constrict in presence of edema or cough</li> <li>• Minimizes possibility of tube dislodgment</li> </ul>
<b>Comfortable and Convenient</b>	<ul style="list-style-type: none"> <li>• Minimizes time required to change trach holder</li> <li>• Minimizes time for tracheostomy care as the fasteners are moisture resistant</li> </ul>	<ul style="list-style-type: none"> <li>• Allows for easy tracheostomy tube care</li> </ul>
<b>Latex-Free</b>	<ul style="list-style-type: none"> <li>• Eliminates patient and care-giver latex reactions</li> </ul>	<ul style="list-style-type: none"> <li>• Prevents latex hypersensitivity</li> </ul>