

# Dale® Hold-n-Place®

## Foley Catheter Holders

Features	Benefits	Clinical Advantages	Hold-n-Place® Adhesive Patch	Hold-n-Place® Leg Band
<b>EXCLUSIVE FOLEY CATHETER LOCKING TABS</b>	<ul style="list-style-type: none"> <li>• Secures all types of Foley catheters</li> <li>• Ensures Foley catheter security</li> <li>• Re-sealable adhesive and green fastener tab can be repositioned several times</li> <li>• Allows surgeon to secure Foley catheter with traction</li> </ul>	<ul style="list-style-type: none"> <li>• No fear of choosing the incorrect size.</li> <li>• Catheter cannot slide out of position.</li> <li>• Reduces risk of urethral erosion.</li> <li>• Reduces the risk of traumatic removal</li> <li>• Stabilizes Foley catheter to prevent pistoning and accidental dislodgement.</li> <li>• Simplifies Foley catheter care</li> <li>• Easy patient monitoring</li> </ul>	<ul style="list-style-type: none"> <li>✓</li> <li>✓</li> <li>✓</li> <li>✓</li> <li>✓</li> <li>✓</li> </ul>	<ul style="list-style-type: none"> <li>✓</li> <li>✓</li> <li>✓</li> <li>✓</li> <li>✓</li> <li>✓</li> </ul>
<b>LOW PROFILE</b>	<ul style="list-style-type: none"> <li>• Designed to reduce the likelihood of the tube being caught on bed or clothing</li> </ul>	<ul style="list-style-type: none"> <li>• Prevents tube dislodgement</li> </ul>	<ul style="list-style-type: none"> <li>✓</li> </ul>	<ul style="list-style-type: none"> <li>✓</li> </ul>
<b>NO HARD PLASTIC PARTS</b>	<ul style="list-style-type: none"> <li>• Nothing to dig into fragile skin</li> </ul>	<ul style="list-style-type: none"> <li>• Reduces likelihood of skin irritation and skin breakdown</li> </ul>	<ul style="list-style-type: none"> <li>✓</li> </ul>	<ul style="list-style-type: none"> <li>✓</li> </ul>
<b>VERSATILE</b>	<ul style="list-style-type: none"> <li>• Universally Sized</li> </ul>	<ul style="list-style-type: none"> <li>• One size fits all Foley catheters</li> <li>• Can be secured at the y-port of the Foley catheter</li> <li>• Can be secured on the drainage bag tubing</li> <li>• One Adhesive SKU fits all patients' leg size</li> <li>• Two Leg Band SKU's fit all patients' leg size</li> <li>• Can be applied to abdomen<sup>†</sup></li> </ul>	<ul style="list-style-type: none"> <li>✓</li> <li>✓</li> <li>✓</li> <li>✓</li> <li>✓</li> <li>✓</li> </ul>	<ul style="list-style-type: none"> <li>✓</li> <li>✓</li> <li>✓</li> <li>✓</li> <li>✓</li> <li>✓<sup>†</sup></li> </ul>

<sup>†</sup>330 Bariatric leg band/waistband



[www.dalemed.com](http://www.dalemed.com)

Dale® & Hold-n-Place® are registered trademarks of Dale Medical Products, Inc. in the USA & EU.  
©2017 Dale Medical Products, Inc.

CA-HNP REV E

# Dale® Hold-n-Place®

Foley Catheter Holders

Features	Benefits	Clinical Advantages	Hold-n-Place® Adhesive Patch	Hold-n-Place® Leg Band
<b>RELIABLE ADHESIVE BASE</b>	<ul style="list-style-type: none"> <li>• Skin-friendly, clear and breathable</li> </ul>	<ul style="list-style-type: none"> <li>• Water-resistant</li> <li>• Writable base surface to record patient data</li> <li>• Leaves no residue on patient's skin</li> <li>• Adherence to skin for up to three days</li> </ul>	<ul style="list-style-type: none"> <li>✓</li> <li>✓</li> <li>✓</li> <li>✓</li> </ul>	
<b>SOFT, SUPPLE STRETCH MATERIAL</b>	<ul style="list-style-type: none"> <li>• Tape-free method ideal for long-term securement</li> </ul>	<ul style="list-style-type: none"> <li>• Does not constrict blood flow of superficial or deep veins</li> <li>• Maintains position on the leg, whether ambulatory or bedridden</li> <li>• Allows leg band to be moved from leg to leg to decrease meatal irritation and skin trauma</li> <li>• Eliminates skin irritation and breakdown from repeated tape application</li> </ul>		<ul style="list-style-type: none"> <li>✓</li> <li>✓</li> <li>✓</li> <li>✓</li> </ul>
<b>IMPROVES PATIENT COMFORT</b>	<ul style="list-style-type: none"> <li>• Allows for easy relocation of holder</li> </ul>	<ul style="list-style-type: none"> <li>• Can be moved/rotated to minimize pressure points</li> <li>• Allows for normal daily activities with no fear of Foley catheter pullout</li> </ul>	<ul style="list-style-type: none"> <li>✓</li> <li>✓</li> </ul>	<ul style="list-style-type: none"> <li>✓</li> <li>✓</li> </ul>
<b>SHORT-TERM SECUREMENT</b>			<ul style="list-style-type: none"> <li>✓</li> </ul>	<ul style="list-style-type: none"> <li>✓</li> </ul>
<b>LONG-TERM SECUREMENT</b>				<ul style="list-style-type: none"> <li>✓</li> </ul>
<b>NOT MADE WITH NATURAL RUBBER LATEX</b>	<ul style="list-style-type: none"> <li>• Eliminates patient and care-giver latex reactions</li> </ul>	<ul style="list-style-type: none"> <li>• Prevents latex hypersensitivity reactions</li> </ul>	<ul style="list-style-type: none"> <li>✓</li> </ul>	<ul style="list-style-type: none"> <li>✓</li> </ul>

US Patent #5,941,856 - Leg Band



[www.dalemed.com](http://www.dalemed.com)

Dale® & Hold-n-Place® are registered trademarks of Dale Medical Products, Inc. in the USA & EU.  
©2017 Dale Medical Products, Inc.

CA-HNP REV E