

SIZING CHART

| Product | Product Number | Description |
|------------------------------------|----------------|---|
| NasoGastric Tube Holders | 160 | Tube Holder, fits most |
| | 161 | Tube Holder, fits small |
| Tracheostomy Tube Holders | 200 | Dale® One Piece 1" wide neckband, fits 9.5" - 18" neck (24-46cm) |
| | 240 | 240 Blue® 1" wide neckband, fits up to 19.5" neck (49cm) |
| | 241 | 241 PediPrints® 1" wide neckband, fits up to an 18" neck (46cm) |
| | 242 | 242 PediPrints® 0.75" wide neckband, fits up to a 9" neck (23cm) |
| | 244 | 244 Blue® 1.5" wide neckband, fits up to 28" neck (71cm) |
| Endotracheal Tube Holders | 270 | Stabilock, fits endotracheal tubes 7.0-10.0mm |
| | 273 | Adhesive Base for Stabilock, one size fits most |
| Foley Catheter Tube Holders | 150 | Hold-n-Place® Adhesive Patch, One size fits most |
| | 316 | Hold-n-Place® Leg Band, fits up to 30" (76cm), One size fits most |
| | 330 | Hold-n-Place® Bariatric Leg Band/Waistband, fits up to 56"(142cm) |
| ACE Connector® | 475 | Access Controller for Enteral feeding, fits most |
| ACE Connector® with ENFit® | 485 | Access Controller for Enteral feeding, fits most ENFit connectors |
| Nasal Dressing Holders | 600 | One size fits most |
| Transducer Holders | 620 | Fits up to three standard transducers, one size fits most |
| | 622 | Fits up to two large transducers, one size fits most |
| Bendable ArmBoards | 650 | Large, 9" x 3.5" (23cm x 12cm) |
| | 651 | Medium, 5.38" x 2.5" (13.7cm x 6.4cm) |
| | 652 | Small, 4.5" x 1.75" (11.4cm x 4.4cm) |
| | 653 | PediPrints® X-Small, 4.25" x 1" (10.8cm x 2.5cm) |
| IV-ARMOR® | 675 | One size fits most |
| Abdominal Binders | 525 | 2 panel, 6" wide stretches to fit 28-50" (72-127cm) |
| | 410 | 3 panel, 9" wide, stretches to fit 30-45" (76-114cm) |
| | 411 | 3 panel, 9" wide, stretches to fit 46-62" (117-157cm) |
| | 418 | 3 panel, 9" wide, stretches to fit 60-75" (152-191cm) |
| | 810 | 4 panel, 12" wide, stretches to fit 30-45" (76-114cm) |
| | 811 | 4 panel, 12" wide, stretches to fit 46-62" (117-157cm) |
| | 818 | 4 panel, 12" wide, stretches to fit 60-75" (152-191cm) |
| | 820 | 4 panel, 12" wide, stretches to fit 72-94" (183-239cm) |
| | 920 | 5 panel, 15" wide, stretches to fit 72-94" (183-239cm) |
| Post-Surgical Bras | 702 | Medium, fits 34-36" B-D (86-92cm) |
| | 703 | Large, fits 36-38" B-D (92-97cm) |
| | 704 | X-Large, fits 38-44" B-D (97-112cm) |
| | 705 | XX-Large, fits 46-54" C-E (117-137cm) |
| | 706 | XXX-Large, fits 54-60" E-I (137-152cm) |

See reverse for additional products

SIZING CHART

| Product | Product Number | Description |
|---|----------------|---|
| Hold-n-Place® Catheter Securement | 848 | IV, Arterial & Mid-Line Securement Kit w/ CHX Transparent Dressing |
| | 849 | IV, Arterial & Mid-Line Catheter Securement Kit w/ Transparent Dressing |
| | 850 | IV, Arterial & Mid-Line Catheter Securement |
| | 852 | CVC, PICC & Arterial Sheath Securement Device |
| | 853 | CVC, PICC & Arterial Sheath Securement Kit w/ Transparent Dressing |
| | 854 | CVC, PICC & Arterial Sheath Securement Kit w/ CHX Transparent Dressing |
| Hold-n-Place® General Purpose Securement | 950 | Small, fits 1.5mm-4.5mm tubing |
| | 951 | Medium, fits 3mm-8mm tubing |
| | 952 | Large, fits 6mm-13mm tubing |