


Dale Hold-n-Place® Foley Catheter Holders secure the catheters to reduce the likelihood of traumatic removal, urethral erosion, bladder spasms and trauma.

Dale[®]
Always Better



Hold-n-Place® Foley Catheter Holders

WHEN IT COMES TO SECURING FOLEY CATHETERS, THERE'S NO REASON TO COMPROMISE.


In addition to saving valuable nursing time, the Dale Hold-n-Place Foley Catheter Holders stabilize the indwelling catheter, reducing the likelihood of traumatic removal, without restricting patient's movement or irritating the skin.

- The leg band, made entirely of soft stretch material, distributes compression evenly around the thigh to minimize circumferential pressure. Ideal for long-term catheterization.†
- The adhesive patch, with a skin-friendly adhesive base, is ideal for short-term catheterized patients.
- The locking devices speed up application and simplify adjustment.
- The holders can be used with equal effectiveness on all Foley catheters including Teflon, silver and silicone coated catheters
- Leg band/waist band ideal for bariatric patients
- Not made with natural rubber latex.

SIMPLIFIED CLOSURE SYSTEM

All Dale Hold-n-Place Foley Catheter Holders feature patented locking systems that securely anchor the catheter in place, reducing the risk of urethral erosion, bladder spasms and trauma. The locking devices speed up application and simplify adjustment.

The holders can be used with equal effectiveness on all Foley catheters including Teflon®, silver and silicone coated catheters.



ORDERING INFORMATION

Product Number	Size	Quantity
150	Adhesive Patch, one size fits most	50 per box
316	Leg band, fits up to 76cm(30")	10 per box
330	Leg band/Waist band, fits up to 142cm(56")	10 per box

†More than 35 million uses over 30 years with no adverse incidences reported to Dale. (Data on file.)

CMS Guidelines: As of October 2008, Medicare will no longer reimburse hospitals for hospital acquired complications such as a catheter-associated urinary tract infections. HAIs are expensive. Make Dale Hold-n-Place Foley Catheter Holders part of your hospital's protocol.

**WOCN recommends stabilizing foley catheters with an adhesive anchoring system or a catheter strap.

***CDC states indwelling catheters should be properly secured after insertion to prevent movement and urethral traction.

***OSHA recommends securing medical catheters.

Dale
www.dalemed.com

Dale® & Hold-n-Place® are registered trademarks of Dale Medical Products, Inc. in the USA & EU.
Teflon® is a registered trademark of E.I. du Pont de Nemours and Company or its affiliates.
©2020 Dale Medical Products, Inc.

PS-HNP REV E